



A Class
Apart

GHK-Cu

50mg pen, reconstituted with 3mL of bacteriostatic water



WHAT IT DOES

- Powerful skin rejuvenation and collagen support
- Supports wound healing and tissue repair
- Promotes hair growth and scalp health
- Antioxidant and anti-ageing peptide

HOW IT IS USED



Taken daily



Subcutaneous injection



Ideal for skin rejuvenation, anti-ageing and recovery



Can be used any time of day



6–12 weeks on, then 4 weeks off

DOSAGE PROTOCOL

Concentration: 50mg / 3mL = 16.67mg/mL

UNITS / CLICKS	VOLUME	DOSE
3 units / clicks	0.10mL	1.67mg
6 units / clicks	0.20mL	3.33mg
9 units / clicks	0.30mL	5mg
12 units / clicks	0.40mL	6.67mg
15 units / clicks	0.50mL	8.33mg
18 units / clicks	0.60mL	10mg

All dose conversions assume the peptide has been reconstituted exactly to 50mg in 3mL.



Skin Rejuvenation / Anti-Ageing:
1–2.33mg (2.33mg optimal)



Hair Growth / Scalp Repair:
1–2.33mg (2.33mg optimal)



Wound Healing / Post-Procedure:
5mg




Duration / Cycle:
6–12 weeks on, then 4 weeks off



Cycle Break Support Option:
Glow Stack

POSSIBLE SIDE EFFECTS

- Mild irritation or redness at injection site
- Occasional nausea, headache or temporary digestive upset



WhatsApp 05262 78691